

VINYL GENTLE TOUCH CLEAR POWDER FREE



VINYL

CLEAR

POWDER FREE

S-XL

AMBIDIEXTROUS

NON- STERILE

High sensitivity glove. It offers protection against specific chemicals risks and against solvents in sporadic contact. It offers a wide variety of general uses highlighting medical use, examination, lab, cleaning, gardening, etc.

Characteristics:

Suitable for latex allergic people.	Reduces fatigue.
Clear, powder free, single use, ambidextrous	High sensitivity
Rolled Edge, protection against spills	For application requiring a continuous change of glove

Packaging:

- 100 pieces par dispenser, with easy opening for a hygienic access to the glove.
- Packaging contains: 10 dispensers of 100 units
- Dimension: 210 x 120 x 60 mm



Labelling: Manufacturing date, expiration date, batch number, reference, and EAN code. Product name in several languages: Spanish, English, French, German, Italian, Portuguese.

Ce marking, Protection Pictograms, Safety Standards and legislation. Size and units, size guide.

Normative:

- Category III. Ce mark according to Regulation (UE) 2016/425
- Class I. Royal Decree 745/2017 which regulates sanitary products, which transposes Directive 93/42 EEC Medical Device.
- Food Normative 10/2011 of 14 January 2011 on plastics materials and articles intended to come into contact with food.
- ISO 16604 Bacteriophage. Clothing for protection against contact with blood and body fluids -- Determination of resistance of protective clothing materials to penetration by blood-borne pathogens -- Test method using Phi-X 174 bacteriophage. ASTM F1671. Resistance to the penetration of blood and bacteriophage virus.
- EN ISO 374-1:2016/A1 :2018. Standard EN 374-2. Point 7. Standard EN374-4, Point 5.
- EN ISO 374-5:2016. Protection against bacteria, fungi and virus.
- EN 16523-2:2015+A1:2018. Resistance to permeation by chemicals.
- EN ISO 21420:2020

RUBBEREX SPAIN S.L

+34 96 166 7124

info@rubberex.es

www.rubberex.es

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España

VINYL GENTLE TOUCH CLEAR POWDER FREE



- EN455 Medical gloves for single use Part1 Requirements and testing for freedom from holes. Part 2 - Requirements and testing for physical properties. Part 3 - Requirements and testing for biological evaluation. Part 4 - Requirements and testing for shelf life determination
- UNE-EN ISO 21171:2007 (Point 5.2. UNE-EN 455-3). Powder free glove, average <2 mg

CHEMICAL RESISTANCE.

PERMEATION. EN ISO 374-1:2016/A1:2018

Gloves are classified in terms of time of passage, according to each individual chemical for which the glove resists permeation.

The gloves are tested against the following chemicals:

CHEMICAL	LEVEL	CODE LETTER
40 PER CENT SODIUM HYDROXIDE	6	K

Performance Level:

Time (min)	> 10	> 30	> 60	> 120	> 240	> 480
Performance Level	1	2	3	4	5	6

DEGRADATION

Degradation % (DR) is determined for each chemical made in the permeation.

(DR) K: 40% Sodium Hydroxide = 30.7%

The Gentle Touch, Vinyl glove PPE should be used to protect the user's hand against chemical risks products mentioned above) and against microbiological risks (bacteria, fungi and virus).

The present PPE should be used against other risks than those described above. The use of this PPE should be evaluated correctly depending on the place of work

RUBBEREX SPAIN S.L

+34 96 166 7124

info@rubberex.es


www.rubberex.es

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España

VINYL GENTLE TOUCH CLEAR POWDER FREE



EN ISO 374-1:2016:A1:2018 / EN ISO 374-5:2016



Gloves Dimensions	S	M	L	XL
Palm Width mm	85 ± 5 mm	95 ± 5 mm	110 ± 5 mm	110 ± 5 mm
Length mm	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm
Palm thickness mm (± 0.02)	0.07	0.07	0.07	0.07
Finger thickness mm (± 0.02)	0.11	0.11	0.11	0.11
Cuff thickness mm (± 0.02)	0.06	0.06	0.06	0.06
Weight gr (± 0.3)	4.2	4.5	4.7	4.9

REF	Size	Size	Size	EAN
0422	Natural	S	7	8435234004227
0423	Natural	M	8	8435234004234
0424	Natural	L	9	8435234004241
0425	Natural	XL	10	8435234004258

Dispenser/ Inner



Logistic:

inner: 100 gloves

Carton: 10 inners per carton

Pallet: 70 cartons

REV 29.jul 23

RUBBEREX SPAIN S.L

+34 96 166 7124

info@rubberex.es

www.rubberex.es

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España