

# LATEX GENTLE TOUCH CLEAR POWDERED



LÁTEX

CLEAR

POWDERED

XS-XL

AMBIDEXTROUS

NON- STERIL

High sensitivity glove. Provides protection against specific chemical risks and solvents in sporadic contact. Latex gloves are suitable for so much purposes as for example: medical use, examination, hygiene and cleaning, hairdresser, tattoo artists, etc.

## Characteristics:

<b>Not suitable for latex allergic people.</b>	<b>Low-Level latex protein.</b>
<b>Clear, powdered, single use, ambidextrous.</b>	<b>High sensitivity and ergonomic.</b>
<b>Rolled Edge, protection against spills</b>	<b>Flexibility and adaptability.</b>

## Packaging:

- 100 pieces par dispenser, with pre-cut centre circle.
- Packaging contains: 10 dispensers of 100 units
- Dimension: 210 x 60 x 120 mm



**Labelling:** Manufacturing date, expiration date, batch number, reference and EAN code. Product name in several languages: Spanish, English, French, German, Italian, Portuguese.

CE marking, protection pictograms, safety standards and legislation. Size, size guide and units.

## Normative:

- Category III. Marca CE mark according to Regulation (EU) 2016/425
- Class I. Royal Decree 2017/745 which regulates sanitary products, which transposes Directive 93/42 CEE.
- Food Normative 10/2011 of 14 January 2011 on plastic material and articles intended to come into contact with food.
- EN ISO 374-1:2016. Standard EN 374-2. Point 7. Standard EN374-4. Point 5.
- EN ISO 374-5:2016
- EN420:2003+A1:2009
- EN455-1-2-3

RUBBEREX SPAIN S.L

+34 96 166 7124

[info@rubberex.es](mailto:info@rubberex.es)

[www.rubberex.es](http://www.rubberex.es)

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España

# LATEX GENTLE TOUCH CLEAR POWDERED



## CHEMICAL RESISTANCE

### PERMEATION. EN ISO 374-1:2016/A1:2018

Gloves are classified in terms of time of passage, according to each individual chemical for which the glove resists permeation:

CHEMICAL	LEVEL	WORD
40% SODIUM HYDROXIDE	2	K

Performance Level:

Time (min)	> 10	> 30	> 60	> 120	> 240	> 480
Performance Level	1	2	3	4	5	6

## DEGRADATION

Degradation % (DR) is determined for each chemical made in the permeation.

(DR) K: 40% Sodium Hydroxide = 4.53%

The “Gentle Touch Latex”, Latex glove PPE should be used to protect the user’s hand against chemical risks (products and levels mentioned above) and against microbiological risks.

The present PPE should never be used against other risks than those described above. The use of his PPE should be evaluated correctly depending on the place of work.

EN ISO 374-1:2016/A1:2018  
Type C

EN ISO 374-5:2016



K



CAT III

RUBBEREX SPAIN S.L

+34 96 166 7124

[info@rubberex.es](mailto:info@rubberex.es)

[www.rubberex.es](http://www.rubberex.es)

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España

# LATEX GENTLE TOUCH CLEAR POWDERED



Glove Dimension	XS	S	M	L	XL
Palm width mm	81 ± 5 mm	85 ± 5 mm	95 ± 5 mm	110 ± 5 mm	110 ± 5 mm
Length mm	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm
Palm thickness mm (± 0.02)	0.09	0.09	0.09	0.09	0.09
Finger thickness mm (± 0.02)	0.10	0.10	0.10	0.10	0.10
Cuff thickness mm (± 0.02)	0.07	0.07	0.07	0.07	0.07
Weight gr. (± 0.3)	4.5	4.8	5.0	5.3	5.5

REF	Colour	Size	Size	EAN
214	Clear	XS	6	8435234502143
194	Clear	S	7	8435234501948
195	Clear	M	8	8435234501955
196	Clear	L	9	8435234501962
197	Clear	XL	10	8435234501979

## Dispenser



## Logistics:

Dispenser: 100 gloves

210 mm x 60 mm x 120 mm

Carton: 10 dispensers par cartón

310 mm x 250 mm x 250 mm

Pallet: 84 export carton

Rev 03 may 2023

RUBBEREX SPAIN S.L

+34 96 166 7124

[info@rubberex.es](mailto:info@rubberex.es)

[www.rubberex.es](http://www.rubberex.es)

Polígono Industrial Masía del Baló. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valencia. España