

GENTLE TOUCH NATURAL LATEX 5.0 GR POWDERED



LATEX

NATURAL

POWDERED

S-L

AMBIDEXTROUS

NON-STERILE

High sensitivity glove. It offers protection against specific chemical hazards and sporadic contact solvents. Its uses include health, examination and examination, hygiene and cleaning, hairdressing, tattoo artists, etc.

Characteristics:

Not suitable for latex allergy sufferers.	Low level of latex protein.
Natural, powdered, single-use, ambidextrous	Is sensitive and ergonomic.
Cuff with flange, spill protection.	Flexibility and Adaptability.

Presentation:

- Dispenser box of 100 units, with petrochelated opening for hygienic access to the glove.
- Packaging contains: 10 boxes of 100 units
- Dimension: 210 x 60 x 120 mm



Labelling: Expiry date, date of manufacture, batch number, reference, barcode. Product name in several languages: Spanish, English, French, German, Portuguese, Italian.

CE marking, Protection pictograms, Safety standards and legislation. Size and units, size guide.

Regulations:

- Category I. CE marked according to Regulation (EU) 2016/425
- Food regulation 10/2011 on plastic materials intended to come into contact with food
- EN420:2003+A1:2009

The PPE type LATEX GLOVE "GENTLE TOUCH" should be used to protect the wearer's hand against chemical hazards (products and levels mentioned above) and against microbiological hazards.

This PPE should never be used against risks other than those described above. The use of this PPE should be evaluated correctly depending on the workplace.

RUBBEREX SPAIN S.L

+34 96 166 7124

info@rubberex.es

www.rubberex.es


Masía del Baló Industrial Estate. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valence. Spain

GENTLE TOUCH NATURAL LATEX 5.0 GR POWDERED



Glove Dimensions	S	M	L	XL
PALM WIDTH	85 ± 5 mm	95 ± 5 mm	110 ± 5 mm	110 ± 5 mm
LONGITUDE	Min 240 mm	Min 240 mm	Min 240 mm	Min 240 mm
Palm thickness mm (± 0.02)	0.09	0.09	0.09	0.09
Finger thickness mm (± 0.02)	0.10	0.10	0.10	0.10
Cuff thickness mm (± 0.02)	0.07	0.07	0.07	0.07
Weight gr. (± 0.3)	4.8	5.0	5.3	5.5

REF	Colour	Size	T. numerical	EAN
194	Natural	S	7	8435234501948
195	Natural	M	8	8435234501955
196	Natural	L	9	8435234501962



Logistics:
 Dispenser: 100 gloves
 210 x 70 x 120 mm
 Box: 10 dispensers per
 box
 360 x 250 x 220 mm
 Pallet: 84 boxes

Rev 10 L 2026

RUBBEREX SPAIN S.L

+34 96 166 7124

info@rubberex.es

www.rubberex.es

Masía del Baló Industrial Estate. C/Mas Vallesa de Mandor N8 Nave12. CP. 46190. Ribarroja del Turia. Valence. Spain